


VWR Scientific



40061-000

Pneumat . . . Air-operated

40061-000 PRESS, Pneumat, BUEHLER (AB No. 20-1361) —

For compression mounting of samples with either thermosetting or thermoplastic materials, in powder or premolded form. The Pneumat Press is air operated from an 80 psi line supply. It provides for outer ejection and "in-position" loading of samples and molding materials. Push-button controls for power, heater and ram simplify operation. There is also a timer and a ram force gauge for reading molding pressure. Panel pressure regulator lets operator alter line supply pressure to suit mold size.

Equipped with built-in air cylinder, rapid mold closure, 1" (2.54 cm.) screw base mold assembly, slip-on cooler and automatic thermostat controlled heater; pair of gloves, samples of Bakelite® and mold release fluid. For 115 volts, 60 Hz.